


# The Clean Cooler Concept™

## PWII-03



The Latest Innovation From



**PURE WATER TECHNOLOGY®**  
DIVISION

# Pure Water II-O<sub>3</sub>

## System Performance Characteristics

- 5-Stage Purification Process
- Microprocessor Control
- Activated Oxygen Injection™
- Self-Sanitizing™
- 50 G.P.D. Thin Film Composite Reverse Osmosis Membrane

### Options:

- Hot Water
- TDS (Total Dissolved Solids) Monitor

## Pure Water II-O<sub>3</sub> vs Bottled Water

- Continuous Supply of Water
- No Heavy Jugs to Lift
- No Storage Issues
- Contaminants are Eliminated in the Purification Process
- Activated Oxygen Injection™
- Self Sanitizing™
- Bottles Must be Changed
- 5-Gallon Bottles Weigh 40 lbs.
- 5-Gallon Bottles Take Up Floor Space
- Bacterial Contaminants can be Introduced When Bottles are Changed
- Bottle Cooler MUST BE Sanitized Regularly



Microprocessor Control

TDS Monitor

Hot, Room Temperature and Cold Faucets

State-of-the-Art Reverse Osmosis and Filtration System

Proudly Distributed By: